




1 2 3 4 5 6 7 8 9

A  
B  
C  
D  
E  
F



NOT TO SCALE

Sheet 1 of 2

NOTE: THIS INFORMATION IS PROPRIETARY. REPRODUCTION OF THIS MATERIAL IS PROHIBITED WITHOUT PRIOR WRITTEN AUTHORIZATION FROM COLE-PARMER.	
Supplier:	<b>COMPONENT</b>
Supplier Part No:	
Cole-Parmer Item:	<b>MASTERFLEX® L/S® Y-TUBING ASSEMBLY</b>

**COMPONENT NAME: MASTERFLEX® L/S® Y-TUBING ASSEMBLY**


Item Key	Part Description	CP Part #	Total Qty in Assembly	Length of Tubing	Notes:
1	Tubing, Masterflex L/S		2	17.0" (+1.0"/-0.0")	
2	Tubing, Masterflex L/S		2	32.0" (+1.0"/-0.0")	
3	Y-Fitting		2		Y Fitting with appropriate banding to ensure engagement of tubing to the fitting.
	Overall the formulations that we would be interested in would be: PharMed BPT and BioPharm Plus				

**Assembly Notes: 1. Attach tubing segments to Y-fitting inserts as shown. 2. Complete with appropriate banding to ensure a leakproof connection when under tubing rated pressures.**



1 2 3 4 5 6 7 8 9

A  
B  
C  
D  
E  
F



NOT TO SCALE

Sheet 1 of 2

NOTE: THIS INFORMATION IS PROPRIETARY. REPRODUCTION OF THIS MATERIAL IS PROHIBITED WITHOUT PRIOR WRITTEN AUTHORIZATION FROM COLE-PARMER.	
Supplier:	<b>COMPONENT</b>
Supplier Part No:	
Cole-Parmer Item:	<b>MASTERFLEX® I/P® Y-TUBING ASSEMBLY</b>

**COMPONENT NAME: MASTERFLEX® I/P® Y-TUBING ASSEMBLY**

Item Key	Part Description	CP Part #	Total Qty in Assembly	Length of Tubing	Notes:
1	Tubing, Masterflex I/P		2	17.0" (+1.0"/-0.0")	
2	Tubing, Masterflex I/P		2	32.0" (+1.0"/-0.0")	
3	Y - Fitting		2		Y Fitting with appropriate banding to ensure engagement of tubing to the fitting.
	Overall the formulations that we would be interested in would be: PharMed BPT and BioPharm Plus				

**Assembly Notes: Assembly Notes: 1. Attach tubing segments to Y-fitting inserts as shown. 2. Complete with appropriate banding to ensure a leakproof connection when under tubing rated pressures.**